



VITA shade taking



VITA Linearguide 3D-MASTER

The new simplicity in shade taking



Product description

With the VITA Linearguide 3D-MASTER you can determine the correct tooth shade swiftly and accurately. The modern design and systematic structure of the VITA Linearguide enable the appropriate 3D-MASTER shade to be found quickly.

The VITA Linearguide 3D-MASTER facilitates the changeover from traditional shade-taking systems to the precise shade determination and shade communication in the VITA SYSTEM 3D-MASTER.

The VITA Linearguide 3D-MASTER contains the same 29 shade samples as the VITA Toothguide 3D-MASTER.

Advantages

- Classical, linear handling
- Self-explanatory, simple use
- Reliable, fast and precise shade determination for satisfied customers
- Modern design as an expression of perfect esthetics and a clear structure

VITA Linearguide products


VITA SYSTEM 3D-MASTER® shades

Description	Content	Prod. no.
VITA Linearguide 3D-MASTER	1 Box	B363

Printed materials

Title	Description	Order no.
VITA Linearguide 3D-MASTER	Instructions	920-02565
VITA Linearguide 3D-MASTER	information sheet	1532

VITA Toothguide 3D-MASTER[®] with BLEACHED SHADE GUIDE



Product description

VITA Toothguide 3D-MASTER[®]

With the VITA Toothguide 3D-MASTER, all natural tooth shades can be determined accurately and systematically.

Advantages

- Determination of the tooth shade in 3 steps
- Unmistakable determination of the tooth shade
- Possibilities of comparison: Yes/No - Better/Worse decision
- Tooth shades precisely defined by science (miscible)
- Can be learned

Product description

BLEACHED SHADE GUIDE

The VITA BLEACHED SHADE GUIDE enables the shades of bleached teeth to be determined accurately and systematically using the same principle as the VITA Toothguide.

VITA Toothguide products


VITA SYSTEM 3D-MASTER® shades

Description	Content	Prod. no.
VITA Toothguide 3D-MASTER	1 pc	B360
VITA BLEACHED SHADE GUIDE, Shade Group 0	1 pc	B2600MS

Printed materials

Title	Description	Order no.
VITA Toothguide 3D-MASTER	Instructions	961
VITA Toothguide 3D-MASTER	Color communication form	914
Determination of tooth shades	Special reproduction	1283
VITA SYSTEM 3D-MASTER	Brochure	1418

VITA Bleachedguide 3D-MASTER®



Product description

The VITA Bleachedguide 3D-MASTER is a tool for the dentist for the planning and control of tooth lightening procedures.

The Composition follows the principles of the VITA SYSTEM 3D-MASTER based on the 3 shade parameters lightness, chroma and hue. In addition to the lightness and chroma levels, however, additional intermediate levels with supplementary shade samples have here been added.

Features/Benefits

- Tooth lightening can be discussed with the patient in realistic terms.
- Uniform color distances which can be clearly distinguished by the human eye
- Intermediate steps are available as shade samples (lightness, chroma)
- Dentist and patient can discuss the desired result on the basis of the shade samples and later check the result obtained.
- Prosthetic restorations can – if necessary – be adapted to the VITA SYSTEM 3D-MASTER.
- At a recall visit it can be checked whether it is necessary to repeat the tooth lightening procedure.

VITA Bleachedguide products

VITA SYSTEM 3D-MASTER® shades

Description	Content	Prod. no.
VITA Bleachedguide 3D-MASTER	1 pc	B361

Printed materials

Title	Description	Order no.
VITA Bleachedguide 3D-MASTER	Instructions	920-26266
VITA Bleachedguide 3D-MASTER	information sheet	1432
VITA Bleachedguide 3D-MASTER	Cardboard (paper) version	1486

VITA classical A1-D4 shade guide



Product description

VITA classical A1 – D4 is a shade guide for determining the shade of teeth.

The arrangement of the shades in the VITA classical family of shades is as follows:

A1 – A4 (reddish-brownish)

B1 – B4 (reddish-yellowish)

C1 – C4 (greyish shades)

D2 – D4 (reddish-grey)

Prod. no. B027C

Tips for shade-taking

- Determine tooth shades if possible in daylight or standardized daylight lamps and not under normal ambient lighting.
- The entire surroundings should be kept free from bright colours. If necessary, ask the patient to remove lipstick, etc., and cover brightly coloured clothing with a grey cape/overall.
- Make your choice quickly, always accept the first decision since the eyes begin to tire after approx. 5–7 seconds.

VITA Easyshade[®] Compact



Product description

VITA Easyshade Compact allows you to determine the correct tooth shade and verify restorations quickly and effectively. The cordless, mobile and lightweight unit is easy to operate and displays shade outputs within a few seconds. The measuring accuracy of the VITA Easyshade Compact spectrophotometer is based on the VITA SYSTEM 3D-MASTER and VITA classical A1 – D4 shade systems to obtain perfect results. With the unique VITA SYSTEM 3D-MASTER all natural tooth shades are systematically determined and completely reproduced.

Features/Benefits

- Fast and reliable shade taking at the push of a button
- Clear vision to the measuring point on the tooth
- Perfectly matched with VITA classical A1-D4 and 3D-MASTER shades
- Independent of lighting conditions and user
- Cordless (battery-powered), mobile and lightweight compact unit
- Quality control thanks to the possibility to verify restorations
- Stores up to 25 shades
- Color communication software for PC
- Extended durability thanks to modern LED technology

VITA Easyshade[®] Compact kit

VITA Easyshade Compact
Prod.no. DEASYCS

Quantity	Content	Material
1	pc	Easyshade Compact Base Unit
1	pc	Easyshade Compact Handpiece
40	pcs	Infection Control Shields
1	pc	CD-Rom (Software)
1	pc	Easyshade Compact Packing
1	pc	mains power supply
1	pc	USB cable
1	pc	Operating Manual
1	pc	Brief instructions for use



Printed materials

Title	Description	Order no.
VITA Easyshade Compact	Instructions for use	1506
VITA Easyshade Compact	Brief instructions for use	1559
VITA Easyshade Compact	Brochure	1533

VITA Easyshade[®] Compact products

Description	Content	Prod. no.
Infection control shield	4 x 40 pcs	D39003
Transport case for Easyshade Compact	1 pc	DEASYCSK

VITA ShadeAssist

Tooth colour determination software



Product description

VITA ShadeAssist is a tooth shade communication software including a database for data of patients, practices and laboratories. The program helps to the user to document and process dental diagnoses of patients and to send them to a dental laboratory or the patient himself.

Features/Benefits

- Software including a database for the communication of tooth shades to be used independent of the shade taking method
- Photo editing and drawing program
- Print or email options
- Bleaching process can be documented
- Tooth contour can be selected independent of the tooth to be prepared

VITA ShadeAssist products

Description	Content	Prod. no.
VITA ShadeAssist	1 CD-Rom	D46800-02