

VITAVM[®]9 kits

VITA SYSTEM 3D-MASTER[®] shades



VITA VM 9 BASIC KIT Basic kit for BASIC layering Prod. no. BV9BK12V2 *

Quantity	Content	Material
26	12 g	BASE DENTINE 1M1-5M3
3	12 g	SUN DENTINE SD1-SD3
3	12 g	CHROMA PLUS CP2, CP3, CP4
2	12 g	ENAMEL ENL, END
1	12 g	NEUTRAL NT
1	12 g	WINDOW WIN
3	12 g	CORRECTIVE COR1, COR2, COR3
1	pc	VITA MODISOL marker
1	50 ml	VITA VM MODELLING LIQUID
1	pk	Firing trays G
1	pk	Fibrous pad firing support
1	pc	Shade indicator
1	pc	VITA Toothguide 3D-MASTER
		Instructions for use

* also available as BASIC KIT SMALL (Prod. no. BV9BKS12V2) in the following shades: 1M1, 1M2, 2M1, 2M2, 2M3, 3L1.5, 3L2.5, 3M1, 3M2, 3M3, 3R1.5, 3R2.5, 4M1, 4M2, 4M3, CP2-CP4, ENL, NT, WIN, COR1, SD1-SD3



VITA VM 9 BUILD UP KIT Build-up kit for BUILD UP layering Prod. no. BV9BUK12 *

Quantity	Content	Material
26	12 g	TRANSPA DENTINE 1M1-5M3
1	50 ml	VITA VM MODELLING LIQUID

* also available as BUILD UP KIT SMALL (Prod. no. BV9BUKS12) in the following 15 shades: 1M1, 1M2, 2M1, 2M2, 2M3, 3L1.5, 3L2.5, 3M1, 3M2, 3M3, 3R1.5, 3R2.5, 4M1, 4M2, 4M3



VITA VM 9 BLEACHED COLOR KIT Very light shades to reproduce bleached teeth Prod. no. BV9BLKV1

Quantity	Content	Material
3	12 g	BASE DENTINE 0M1, 0M2, 0M3
3	12 g	TRANSPA DENTINE 0M1, 0M2, 0M3
1	12 g	ENAMEL LIGHT ENL
1	12 g	NEUTRAL NT
1	12 g	WINDOW WIN
1	50 ml	VITA VM MODELLING LIQUID
1	pc	BLEACHED SHADE GUIDE shade group OM
		Instructions for use

VITA VM 9 ONE COLOR KIT 2M2

One-shade kit for initial tests

Prod. no. BV9OCK2M2V2

Quantity	Content	Material
1	12 g	BASE DENTINE 2M2
1	12 g	ENAMEL ENL
1	12 g	WINDOW WIN
1	5 g	AKZENT glaze Akz25
1	50 ml	VITA VM MODELLING LIQUID
1	20 ml	AKZENT fluid
1	pc	Shade indicator stick
		Instructions for use

* also available with 3M2 (Prod. no. BV9OCK3M2V1)

VITA VM 9 PROFESSIONAL KIT

Additional materials for achieving special shade effects and individualizing the restoration

Prod. no. BV9PKV1 * / **

Quantity	Content	Material
11	12 g	EFFECT CHROMA EC1-EC11
3	12 g	MAMELON MM1-MM3
3	12 g	EFFECT PEARL EP1-EP3
5	12 g	EFFECT OPAL EO1-EO5
11	12 g	EFFECT ENAMEL EE1-EE11
6	12 g	EFFECT LINER EL1-EL6
1	pc	Shade indicator for EFFECT ENAMEL
1	pc	Shade indicator for EFFECT OPAL
1	pc	Shade indicator for EFFECT CHROMA
1	pc	Shade indicator for EFFECT LINER

* also available as PROFESSIONAL KIT SMALL (BV9PKS) in the following 15 shades: EC1, EC4, EC6, EC8, EC9, MM2, EP1, EO2, EE1, EE3, EE7, EE8, EE9, EE10, EE11

** compatible with VITA SYSTEM 3D-MASTER and VITA classical A1-D4 shades



VITA VM 9 GINGIVA KIT

Natural-looking gingiva pastes

Prod. no. BV9GK *

Quantity	Content	Material
5	12 g	GINGIVA G1-G5
1	pc	Shade indicator for GINGIVA

* compatible with VITA SYSTEM 3D-MASTER and VITA classical A1-D4 shades





VITA VM 9 MARGIN KIT

Only for minor corrections in the area of margins

Prod. no. BV9MK *

Quantity	Content	Material
6	12 g	MARGIN M1, M4, M5, M7, M8, M9
1	pc	Shade indicator for MARGIN

* compatible with VITA SYSTEM 3D-MASTER and VITA classical A1-D4 shades



VITA VM 9 ESTHETIC KIT for VITABLOCS for CEREC

Additional materials for customizing VITABLOCS

Prod. no. BV9EKC

Quantity	Content	Material
1		Sample set of VITABLOCS 3D-Master
1	12 g	WINDOW WIN
1	12 g	NEUTRAL NT
2	12 g	ENAMEL ENL, END
1	12 g	EFFECT PEARL EP1
2	12 g	EFFECT ENAMEL EE1, EE10
1	12 g	CORRECTIVE COR1
1	5 g	AKZENT finishing agent Akz 25
1	7 g	Shading paste glaze SP15
1	12 g	EFFECT OPAL EO2
2	12 g	EFFECT CHROMA EC1, EC4
1	12 g	MAMELON MM2
1	50 ml	VITA VM MODELLING LIQUID
1	15 ml	VITA Shading paste liquid
1	pk	Firing trays
1	pk	Platinum pins, 2 pcs
		Instructions for use



VITA VM 9 ADD-ON KIT

Set including individualization materials and accessories

Prod. no. BV9AOK *

Quantity	Content	Material
8	12 g	ADD-ON, ADD1-ADD8
1	12 g	Glaze LT
1	50 ml	VITA VM MODELLING LIQUID
1	20 ml	AKZENT Fluid
1	pk	Fibrous pad firing support
1	pc	Shade indicator for ADD-ON
1	pc	No. 3/0 a & e brush
1		Instructions for use VITA PM 9

* compatible with VITA SYSTEM 3D-MASTER and VITA classical A1-D4 shades